



**Translate**  
**7T research power**  
**into clinical care**  
**MAGNETOM Terra**



# **MAGNETOM Terra – Translate 7T research power into clinical care**

MAGNETOM Terra is designed to let you explore new territories in MRI by enabling powerful 7T research and enhancing clinical care. Uncover a whole new world of clinical insights with double SNR<sup>1</sup> for more precision. Our advanced ultra-high-field (UHF) technology will keep you at the cutting edge of MRI, to attract the brightest minds to your facility, sharpen your competitive edge and strengthen your reputation. It delivers a fertile platform for unlocking research capabilities, publishing new insights first, and setting the pace in diagnostic imaging. Welcome to an exclusive research community. Welcome to a whole new world in MRI.

**Welcome to clinical 7T.**

## **Content**

MAGNETOM Terra at a glance	04
<b>Uncover a whole new world of clinical insights</b>	<b>07</b>
with double SNR for more precision	08
<b>Change the game in UHF business</b>	<b>31</b>
with Siemens' 50% lighter 7T magnet	32
<b>Unlock research beyond clinical limits</b>	<b>37</b>
with 16-channel parallel transmit	38
<b>Join the largest research community</b>	<b>55</b>
with over 75% of all UHF users	56
Service and exchange	60
Technical specifications	64
Siemens Healthineers	66

# MAGNETOM Terra – Translate 7T research power into clinical care

## Unique Dual Mode functionality<sup>2</sup>

- Uncover a whole new world of clinical knowledge with the flexibility to get more from your scanner
- Secure switch between research and clinical operation in less than 7 minutes
- Operating with 2 separate databases for clear research and clinical distinction<sup>2</sup>

## 50% lighter 7T magnet technology<sup>5</sup>

- Released for clinical use in Europe and the U.S.
- Lower weight and cold-shipment for easier integration in clinical environments<sup>5</sup>
- Reduced operating costs thanks to Zero Helium boil-off<sup>6</sup>

## 80/200 gradients, 16-channel pTX<sup>2</sup> and up to 64 channels<sup>3</sup>

- More power for greater diffusion MRI and functional MRI with 80/200 gradients
- Higher homogeneity for challenging body regions with up to 16-channel parallel transmit (in research mode)<sup>2</sup>
- Higher acceleration factors with 64 receive channels<sup>4</sup>

## Double SNR for more precision<sup>1</sup>

- 0.2 mm in-plane resolution to visualize previously unseen structures<sup>7</sup>
- 0.14 cm<sup>3</sup> voxel sizes for metabolic brain mapping<sup>8</sup> (in research mode)<sup>2</sup>
- Submillimeter BOLD fMRI precision to visualize sub-cortical activations<sup>9</sup>



## Explore metabolism

- $^{23}\text{Na}$  MRI and  $^{31}\text{P}$  MRS for metabolic insights in clinical mode<sup>10</sup>
- Dedicated  $^{23}\text{Na}$  head coil and  $^{31}\text{P}$  loop coil
- Broadband RFPA enabling X-nuclei MR with up to 10 nuclei in research mode<sup>2</sup>

## > 75% global market share in 7T technology

- Over 75% of 7T and 100% of vendor-integrated > 7T MRI human scanners worldwide from Siemens
- 7 of 11 leading U.S. hospitals with a 7T, (2018–2019), trust Siemens when they decide for 7T investment<sup>11</sup>
- 73% of ISMRM UHF abstracts in 2018 were based on data from Siemens UHF systems<sup>12</sup>

## The world's largest UHF community

- Largest installed base for exchanging ideas in a strong collaborative network
- An opportunity to enhance your reputation and competitiveness
- Incentive for the brightest minds in the MRI community to work with you



## Uncover a whole new world of clinical insights

Discovering new ground in MRI can help you significantly enhance patient outcomes. Imaging at 7T offers more than double the SNR of 3T. This delivers potential for better lesion conspicuity, faster image acquisition to reduce motion artefacts, and earlier disease detection at submillimeter resolution. MAGNETOM Terra is released for clinical use within Europe and the US. Its Dual Mode lets you switch between clinical and research tasks, unlocking new opportunities and providing a solid, well-founded platform for innovative results.

*“Based on higher resolution, 7T provides new insights into gray and white matter disease in the brain, such as multiple sclerosis, focal cortical dysplasia, and hippocampal sclerosis. Furthermore, functional MR benefits from 7T based on a clinically relevant increase in functional sensitivity and specificity. In musculoskeletal imaging, 7T enhances the visualization of small joint structures and subtle pathologies, such as small meniscal tears, triangular fibrocartilage lesions, and early stages of cartilage degeneration.”<sup>13</sup>*

Professor Siegfried Trattnig  
Director of the MR Centre of Excellence,  
MedUni Wien, Vienna, Austria

# Uncover a whole new world of clinical insights – Double SNR<sup>1</sup> for more precision with clinical applications in Dual Mode

## Dual Mode flexibility<sup>2</sup>

MAGNETOM Terra is the first 7T scanner released for clinical use within Europe and the US. With the release for selected neurological and musculoskeletal scan protocols, it has potential to uncover a whole new world of clinical care. Its unique Dual Mode functionality lets you switch between research and clinical operation, giving you flexibility to get more from your scanner.

## Ultra-fine anatomical resolution

In brain and musculoskeletal MRI, 7T reveals details previously unseen at lower field strengths. For example, cerebral cortex imaging at 0.2 mm in-plane resolution<sup>7</sup> may detect changes in cortical structure indicating early dementia. It also helps visualize cortical microinfarcts and plaques in MS patients and delivers insight into the plaque-vessel relationship and iron accumulation.

## Submillimeter fMRI<sup>9</sup>

The BOLD contrast increases linearly with field strength. In clinical use, this could mean higher precision in oncology compared to 3T applications, for example, through smaller voxel sizes. Potentially, this can increase the accuracy of neurological pre-surgical evaluation of eloquent areas before tumor removal, while keeping scanning times viable<sup>14</sup>.

## Powerful image reconstruction

MAGNETOM Terra delivers improvements in workflow for easier operation and better patient handling. Leveraging the latest syngo MR E12 software platform, it lets you work in the same way as you do with cutting-edge 3T technology. What's more, it comes with the most powerful MaRS (Measurement and Reconstruction System) computing technology ever built.<sup>5</sup>

## Physiology is at your fingertips

MAGNETOM Terra is the first 7T MRI scanner to unleash the full potential of the increased MR signal with multinuclear imaging and spectroscopy in clinical settings. The multinuclear option allows the use of two dedicated coils – a <sup>23</sup>Na head coil and a <sup>31</sup>P loop coil, to explore metabolic insights<sup>10</sup>.





Switch between  
research and  
clinical tasks with  
**Dual Mode<sup>2</sup>**

**Ultra-fine**  
resolution  
to visualize  
details previously  
unseen

Explore  
physiology with  
**multinuclear  
MR<sup>10</sup>**

Ultra-fast image  
reconstruction<sup>5</sup> and  
**syngo MR E12**


**Dual Mode<sup>2</sup>** offers the flexibility to switch from research to clinical tasks

### Clinical Mode

- 1 transmit channel
- 11 kW RF power
- 2  $^1\text{H}$  coils (Head 32, Knee 28)
- 2 MNO coils ( $^{23}\text{Na}$  Head 32,  $^{31}\text{P}$  Flex Loop)
- Neuro and MSK optimized clinical applications


### Research Mode

- single channel and up to 16-channel parallel transmit
- 16 x 2 kW RF power
- Wider range of RF coils
- Whole-Body WIP Applications
- Broadband RFPA for X-nuclei MR (10 nuclei)





- syngo MR E12 software line
- XR Gradients 80/200
- Up to 64 receive channels
- Latest MaRS computer
- 3<sup>rd</sup> order shims




# Clinical Mode – Multiple Sclerosis

Hyperintense MS lesion with hypointense center

**DarkFluid TSE**


0.3 x 0.3 x 3 mm<sup>3</sup>,  
5:59 min



**SWI minIP/phase**

0.2 x 0.2 x 1.2 mm<sup>3</sup>,  
5:38 min


Typical central vein and perivenular demyelination is visible.



Multiple Sclerosis  
with low lesion load



Multiple Sclerosis  
with high lesion load



**Dark Fluid TSE**


0.3 x 0.3 x 3 mm<sup>3</sup>, 5:20 min

**T2 TSE**

0.2 x 0.2 x 3 mm<sup>3</sup>, 5:33 min

**SWI**

0.2 x 0.2 x 1.2 mm<sup>3</sup>,  
5:38 min





# Clinical Mode – Tumor


## 3D SWI of Glioblastoma

3D SWI minIP provides superior assessment of the microvasculature.

*Erwin L. Hahn Institute for MRI, Essen, Germany*


### 3 Tesla

0.85 x 0.72 x 2 mm<sup>3</sup>



### 7 Tesla


0.25 x 0.25 x 1 mm<sup>3</sup>




## **Glioblastoma**

Higher SNR for ultra-high 0.2 mm in-plane resolution for imaging tumor vascularization.

*DKFZ, Heidelberg, Germany*



## Clinical Mode – Stroke




**DarkFluid TSE**


0.3 x 0.3 x 3 mm<sup>3</sup>,  
5:59 min

Small focus of strong cortical/  
subcortical signal abnormality  
in the right precentral gyrus.

*FAU, Erlangen, Germany*




**RESOLVE** 1 x 1 x 3 mm<sup>3</sup>, 1:46 min



**PD FS TSE** 0.2 x 0.2 x 3 mm<sup>3</sup>, 5:14 min



## Clinical Mode – Tumor



**DarkFluid TSE**


0.4 x 0.4 x 3 mm<sup>3</sup>,  
4:22 min

High resolution standard protocols  
for detailed visualisation of  
pathologies, increased tissue  
contrast and high resolution at 7T.

*FAU, Erlangen, Germany*




**SWI minIP** 0.2 x 0.2 x 3 mm<sup>3</sup>, 5:38 min



**T2 TSE** 0.2 x 0.2 x 3 mm<sup>3</sup>, 5:33 min

## Clinical Mode – Healthy volunteer



### Hippocampus imaging

High-resolution imaging of the hippocampus at 0.25 mm in-plane resolution.

*Scannexus, Maastricht, Netherlands*

### T2\* weighted imaging


High-resolution imaging of the brainstem at 0.3 mm in-plane resolution.

*MGH, Boston, USA*






## Clinical Mode – Diffuse axonal injury



3 Tesla PD TSE,  $0.7 \times 0.4 \times 5 \text{ mm}^3$ , TA 2:38 min



3 Tesla T2 TSE,  $0.7 \times 0.4 \times 5 \text{ mm}^3$ , TA 2:38 min



7 Tesla PD PD TSE,  $0.2 \times 0.5 \times 3 \text{ mm}^3$ , TA 3:09 min




7 Tesla T2 TSE,  $0.2 \times 0.5 \times 3 \text{ mm}^3$ , TA 3:09 min

### Diffuse axonal injury

The higher sensitivity at 7T reveals hemosiderin from traumatic brain injury in PD images.

*Erwin L. Hahn Institute for MRI, Essen, Germany*

## Clinical Mode – Healthy volunteer



Clear identification of anatomical structures with increased tissue contrast and high resolution at 7T.

0.6 x 0.6 x 0.6 mm<sup>3</sup>, 13:45 min

FAU, Erlangen, Germany



# Clinical Mode – Stroke

Time of Flight (ToF) with 400 micron isotropic resolution reveals smallest vessels in the brain. The higher the signal and the longer the T1 at 7T are, the higher the quality of the Maximum Intensity Projection (MIP) gets.

Coronal



Sagittal



Axial

Visualize smallest vessels with  
0.4 mm isotropic resolution.

0.4 x 0.4 x 0.4 mm<sup>3</sup>, 8:09 min

FAU, Erlangen, Germany







## Clinical Mode – Healthy volunteer

### 0.2 mm in plane resolution

The basal ganglia Caudate, Putamen and Globus Pallidus can be differentiated.  
Enlarged sections: cortical veins can be depicted.



SWI 0.2 x 0.2 x 1 mm<sup>3</sup>, 10:59 min

# Clinical Mode – Enchondroma

Fine structure visible in the lesion with different contrasts.



PD TSE FS  $0.2 \times 0.2 \times 2.5 \text{ mm}^3$ , 3:15 min



T2 TSE  $0.3 \times 0.3 \times 2 \text{ mm}^3$ , 3:24 min



3D DESS  $0.5 \times 0.5 \times 0.5 \text{ mm}^3$ , 3:43 min



## Clinical Mode – Healthy volunteer

Clear delimitation of anatomical structures, such as ligaments, vessels or cartilage.

### T1 SE

0.2 x 0.2 x 2.5 mm<sup>3</sup>,  
4:05 min



### T1 qSE

0.3 x 0.3 x 2.5 mm<sup>3</sup>,  
7:21 min





## Clinical Mode – Healthy volunteer



**PD qTSE FS**

0.2 x 0.2 x 2.5 mm<sup>3</sup>,  
3:15 min



**T1 FL3D WE**


0.5 x 0.5 x 0.5 mm<sup>3</sup>,  
4:35 min

## Clinical Mode – Healthy volunteer

High resolution fiber tracking  
with SMS RESOLVE at 7T  
and syngo.via Frontier<sup>15</sup>

1.4 mm isotropic, 29:22 min


FAU, Erlangen, Germany




CAIPIRINHA acceleration for  
the SPACE pulse sequence at 7T,  
acceleration: 3x2

0.7 mm isotropic, 6:52 min

FAU, Erlangen, Germany




# Clinical Mode – Healthy volunteer



<sup>31</sup>P FID of the human calf  
with NOE enhancement

Non selective 1:42 min

FAU, Erlangen, Germany



<sup>31</sup>P CSI of the human liver

15 ml voxel volume, 4:42 min


FAU, Erlangen, Germany

## Clinical Mode – Healthy volunteer

Additional metabolic  
information with  $^{23}\text{Na}$  UTE,  
fused with T2 TSE

3 mm isotropic, 7:38 min

FAU, Erlangen, Germany






## Clinical Mode – Tumor

Additional metabolic  
information with  $^{23}\text{Na}$  UTE,  
fused with T2 TSE

3 mm isotropic, 7:38 min

FAU, Erlangen, Germany





## Change the game in UHF business

Medical research funding has stagnated in the last decade. Ensuring your high-end MRI endeavors have the right business impact is crucial in today's competitive environment. MAGNETOM Terra is the result of over 25 years of Siemens UHF innovations, culminating in the design of a brand-new, volume-produced 7T magnet. The magnet is 50% lighter than previous generations and supports easier integration into clinical environments. MAGNETOM Terra can help you become more competitive, while making a tangible difference to clinical care, research – and your business.

*"When you talk to other people in the field, it is clear that Siemens has by far the greatest expertise in ultra-high-field imaging."<sup>13</sup>*

Professor Rainer Goebel  
University of Maastricht & scannexus,  
Maastricht, The Netherlands

# Change the game in UHF business with Siemens Healthineers' 50% lighter 7T magnet<sup>5</sup>

## **Innovative magnet technology**

Siemens Healthineers' 7T magnet is a milestone in MR magnet technology. Its unique design and thermally balanced materials minimize physical interactions between core components. The result is 50% lighter than previous generations, with a higher structural stability and a greater fundamental stress capacity. In addition, excellent homogeneity makes for enhanced image quality.

## **Easy clinical integration**

Thanks to the lighter magnet, the scanner can be shipped cold via airfreight. What's more, you benefit from up to 50% faster installation time and ramp-up. Zero Helium boil-off<sup>6</sup> translates into lower lifecycle costs and an improved eco-footprint. All this has the potential to enhance performance, lower resource consumption, improve sustainability, and reduce operating costs.

## **Increased competitiveness**

MAGNETOM Terra<sup>8</sup> can help you broaden research funding opportunities, making your institution stand out as a leader in life sciences. By being at the cutting edge of clinical care and research, you have the opportunity to increase competitiveness for grants, benefit from reduced complexity in clinical trials, and open up potential for clinical imaging reimbursements.

## **Forward-looking technology**

An investment in MAGNETOM Terra is an investment in the future. Siemens is committed to serving the ultra-high-field community – today and tomorrow – with a host of outstanding innovations. From development and production, to service – all of MAGNETOM Terra's key components are delivered from a single reliable partner you can trust, for maximum peace of mind.





Released for  
**clinical use**  
in Europe and  
the U.S.

**50%**  
lighter magnet  
technology<sup>5</sup>

Lower weight and  
**cold-shipment**  
for easy  
integration

**Zero**  
Helium boil-off<sup>6</sup>

## Proven innovations in the development and production of magnet technology

**1980**

The world's first superconducting whole-body MRI

**1994**

The world's first open MRI magnet

**1989**

The world's first 1.5T active-shielded magnet


**1997**

The world's first 3T active-shielded magnet

### Award-winning development and production

Siemens Magnet Technology in Oxford, UK, has received seven Best Factory Awards and seven Queen's Awards for Enterprise in multiple categories, including for processes and design. The facility deploys leading-edge supply chain management methods, and prides itself on reliable, robust production and the highest standards of quality.





#### April 2015 – Installation of Siemens' first 7T magnet

During the 30 years that we have been producing 1.5T and 3T magnets, we have gained extensive engineering skills and well-founded process expertise. This knowledge and experience has led to the development and production of our own 7T magnet.



Universitätsklinikum Erlangen, Germany

*"We are extremely proud at Siemens Magnet Technology to have developed the 7T magnet at the heart of the MAGNETOM Terra. Once again our expert design and process teams have demonstrated how their innovative thinking has led to a product that has pushed forward the boundaries of magnet technology. It is wonderful to see how seamlessly the manufacture of this flagship product has already been integrated into our award winning facility."*

**Ralph Seidler**  
Managing Director,  
Siemens Magnet Technology





## Unlock research beyond clinical limits

When it comes to research, the freedom to push the boundaries is imperative for gaining a competitive edge. For neuroscience and clinical research applications, MAGNETOM Terra delivers ultra-high SNR and up to 16-channel pTX<sup>2</sup> for imaging challenging body regions. In addition, it has up to 64 receive channels<sup>3</sup> for higher acceleration factors<sup>4</sup> and 80/200 gradients for maximum flexibility, as well as a broadband RFPA for the multinuclear option. Moreover, this powerful, reliable scanner supports basic research by helping you develop groundbreaking technologies, set new trends and translate your outcomes into clinical routine.

*“The increased spatial resolution offered by 7T MRI enables us to study fine-grained activation patterns within cortical areas and investigate detailed functional topography of the cerebral cortex in individual human subjects. This will provide us with a deeper understanding of the human brain and its connectomics in healthy and diseased populations.”<sup>13</sup>*

Professor Kamil Ugurbil  
Director of the Center for Magnetic Resonance Research (CMRR),  
Minneapolis, Minnesota, USA

# Unlock research beyond clinical limits with 16-channel parallel transmit<sup>2</sup>

## Enhanced images with pTX<sup>2</sup>

Image quality and speed are key, but inhomogeneities may present challenges, for example, in body MRI. MAGNETOM Terra's pTX technology with up to 16 channels helps you overcome these issues and generate images of excellent quality. This particularly promising technology has the potential to support your own hardware developments.

## More power for your research

MAGNETOM Terra offers a host of cutting-edge research functionalities, providing access to works in progress packages<sup>2</sup> and powerful hardware configurations. 80/200 gradients and up to 64 receive channels deliver enhanced capabilities for your studies. What's more, the scanner gives you the freedom to explore and develop new clinical applications only possible at ultra-high-field strengths.

## Ultra-high resolution spectroscopy<sup>2</sup>


Proton magnetic resonance spectroscopy at 7T not only delivers metabolic information, but also gives accurate anatomical insight. Ultra-high 0.14 cm<sup>3</sup> resolution<sup>8</sup> has the potential to reveal valuable new diagnostic information for clinical applications – including patients with tumors, epilepsy, multiple sclerosis or other neurodegenerative diseases.

## Open platform architecture<sup>2</sup>

MAGNETOM Terra provides a flexible, fertile ground for your own UHF hardware and software developments. For example, Siemens collaboration partners benefit from technical support and direct access to the sequence, the Image Calculation Environment (ICE), and imaging protocols.

## Discover a world beyond anatomy

The broadband RFPA enables the acquisition of spectra and images with up to 10 nuclei<sup>2</sup>. While <sup>23</sup>Na and <sup>35</sup>Cl imaging can give insights to salt balance, <sup>31</sup>P spectroscopy can shed a light on energy metabolism<sup>10</sup>. Use this additional layer of information to dive into a world of metabolic imaging with MAGNETOM Terra.



**16-channel**  
pTX for higher  
homogeneity<sup>2</sup>

**80 / 200**  
gradients and  
64 receive channels<sup>3</sup>  
for more research  
power

**Open**  
**platform**  
architecture  
for own  
developments<sup>2</sup>

**Dive into**  
**Physiology with**  
**X-nuclei MR**  
for up to  
10 nuclei<sup>2,10</sup>






### Transparent fibers

Tracks calculated with spherical deconvolution based on diffusion-weighted EPI acquisitions with 1 mm isotropic resolution covering the whole brain. The high SNR provided by 7T allows resolving crossing fibers in many brain sub-regions.

*Max Planck Institute, Leipzig, Germany*





#### In-vivo histology

*syngo.via* Frontier<sup>15</sup>, the research extension of *syngo.via*, helps bridge the gap in post-processing translational research. Cinematic rendered images based on MR data sets may be used for patient counseling, surgery planning, or teaching purposes.<sup>16</sup>




### Transparent fibers

Tracks calculated with spherical deconvolution based on diffusion-weighted EPI acquisitions with 1 mm isotropic resolution covering the whole brain. The high SNR provided by 7T allows resolving crossing fibers in many brain sub-regions.


*Max Planck Institute, Leipzig, Germany*






#### **In-vivo histology**

*syngo.via* Frontier<sup>15</sup>, the research extension of *syngo.via*, helps bridge the gap in post-processing translational research. Cinematic rendered images based on MR data sets may be used for patient counseling, surgery planning, or teaching purposes.<sup>16</sup>



3 Tesla



7 Tesla



## Functional MRI

Contrast-to-noise ratio maps in resting state fMRI.


*Consortium The Human Connectome Project.*

*CMRR, Minnesota, USA*

*Washington University St. Louis, USA*

*Oxford University, UK*






### Depiction of cortical layers

Post-processed high-resolution anatomical MR data reveals reconstructed surfaces at different cortical depth levels. The inner red surface runs along the white/gray matter boundary. The outer blue surface runs along the outer (pial) boundary of the cortex.


*Scannexus, Maastricht, The Netherlands*



### Tissue segmentation

Delineation of subcortical nuclei in the thalamus and brainstem at 1 mm isotropic resolution. White matter nulled MPAGE (top TA 8:52 min) and gray matter nulled MPAGE (bottom TA 10:38 min).

MGH, Boston, USA




#### High-resolution fMRI

Cortical-layer-specific activation with fMRI at 1 mm isotropic resolution, inflated view. The fMRI visual stimulus was designed to activate a pattern in the shape of the number “7” using the known retinotopic mapping in the human visual cortex.

MGH, Boston, USA

*“Siemens provides the best possible open hardware and software environment to explore these new transmit and receive concepts, all of which have proved essential to allow ultra-high fields to fulfill their potential for the benefit of human health.”<sup>13</sup>*

Professor Lawrence L. Wald, Director  
MGH NMR Core at Martinos Center, Department of Radiology,  
Boston, Massachusetts, USA



### Spine imaging

Ultra-high resolution of the cervical spine using a custom-built spine coil.

*MGH, Boston, USA*



**3 Tesla** 0.5 mm in-plane




**7 Tesla** 0.5 mm in-plane



**7 Tesla** 0.3 mm in-plane






### Body imaging

Left image, 3D VIBE FatSat, right image, thin MIP from the 3D VIBE FatSat. Images acquired using pTX and custom-built coils.

*Erwin L. Hahn Institute for MRI, Essen, Germany*




### Body imaging

Image left, rectum carcinoma imaged with a ce-FLASH ( $0.3 \times 0.6 \times 2 \text{ mm}^3$ , TA 2:14 min). Images acquired using pTX and custom-built coils. Image right, abdominal small bowel imaging with TrueFISP ( $0.4 \times 0.8 \times 2 \text{ mm}^3$ , TA 26 s).

*Erwin L. Hahn Institute for MRI, Essen, Germany*




**Short axis view** FLASH cine retro  
(1 x 1 x 4 mm<sup>3</sup>, GRAPPA 2, TA 18 s).



**Right ventricular output tract** FLASH cine retro  
(1.3 x 1 x 4 mm<sup>3</sup>, SENSE 2, TA 13 s).



**Four chamber view** FLASH cine retro  
(1 x 1 x 4 mm<sup>3</sup>, GRAPPA 2, TA 16 s).





**Two chamber view** FLASH cine retro  
(1.2 x 1 x 4 mm<sup>3</sup>, SENSE 2, TA 13 s).

### Cardiac imaging

Accelerated T1-weighted FLASH acquisitions using a customer-built coil.

*Berlin Ultrahigh Field Facility, Berlin, Germany*



### Hip cartilage transplant


Cartilage transplant visible in 3D DESS ( $0.37 \times 0.74 \times 0.74 \text{ mm}^3$ , TA 5:11), 3D VIBE SPAIR ( $0.19 \times 0.39 \times 0.8 \text{ mm}^3$ , TA 5:58) and T2 MapIt ( $0.25 \times 0.5 \times 2.5 \text{ mm}^3$ , TA 4:47) 14 months after Autologous Chondrocyte Transplantation (ACT) using a customer-built coil.

*Erwin L. Hahn Institute for MRI, Essen, Germany*







Morphologic PD TSE F ( $0.4 \times 0.4 \times 2 \text{ mm}^3$ )



T2 map ( $0.6 \times 0.6 \times 1 \text{ mm}^3$ ), T2 in ms  
More water, disturbed collagen architecture visible



gagCEST image ( $0.7 \times 0.7 \times 3 \text{ mm}^3$ )  
gagCEST asymmetries in [%] lower values, less PG content




Sodium image  
Sodium SNR lower values, less PG content

### Biochemical imaging using CEST


Male patient nine years after Autologous Osteochondral Transplantation (AOT) in the medial femoral condyle.

MedUni Wien, Vienna, Austria





Pre-contrast



Post-contrast



Subtraction

### Breast imaging

High-resolution bilateral breast imaging in short acquisition time. Pre- and post-contrast 3D FLASH SPAIR, 1.4 mm isotropic resolution, TA 1:48 min/series using a customer-built coil.

*NYU Langone Medical Center, New York, USA*



## Join the largest research community

Your reputation plays a pivotal role in your institution's success. MAGNETOM Terra has the power to let you go deeper than ever before, making your research and patient outcomes stand out from the rest. What's more, this leading-edge technology can help you attract the brightest minds to your facility, further enhancing your capabilities. MAGNETOM Terra has the potential to put your organization firmly on the map, offering access to an exclusive network of expertise and broad scope for collaboration and exchange.

*"When we were in a position to order a 7T system, Siemens was the logical choice."<sup>13</sup>*

Professor Peter Jezzard  
Professor of Neuroimaging,  
University of Oxford, Oxford, UK

## Join the largest research community with over 75% of all UHF users

### Enhance your reputation


MAGNETOM Terra helps you achieve your research goals, giving you the opportunity to publish first and become a true opinion leader. This advanced technology has the potential to strengthen your position by attracting the brightest brains to your facility. The scanner lets you deliver previously unseen insights that could improve patient outcomes and further enhance your reputation.

### Expand your network

Even if you are taking your first steps in ultra-high-field imaging, you will never be alone. Siemens has proven expertise in UHF MRI and cultivates links with an extensive network of users. As a result, you benefit from the experience of others and can share your own ideas. Siemens is the global leader in 7T – with a market share of over 75% and more than 25 years of experience in this field.


### Exchange your ideas with peers

When you become part of the Siemens UHF community, you join an exclusive group of outstanding MRI experts. Through collaboration and exchange with other leaders in your field, you can extend your own knowledge and gain deeper insights. Siemens' regular user meetings and an online discussion board are the ideal platforms to interact with your peers.



**73% of ISMRM UHF abstracts** in 2018 were based on data from Siemens UHF systems.<sup>12</sup>






Over **75%**  
of 7T scanners  
deployed  
worldwide are  
from Siemens

Attract and  
retain the  
brightest minds and  
**publish  
first**

**Strong  
network**  
for collaboration  
and peer-to-peer  
exchange



- 7T MRI scanners
- 9.4T, 10.5T and 11.7T MRI scanners

#### UHF systems installed and projects in progress

- 1 Athinoula A. Martinos Center for Biomedical Imaging of MGH, Boston, Massachusetts, USA
- 2 Leibniz Institute for Neurobiology (LIN), Magdeburg, Germany
- 3 Bernard and Irene Schwartz Center for Biomedical Imaging (CBI) of New York University Langone Medical Center, New York City, New York, USA
- 4 Center for MR Research (CMRR), University of Minnesota, Minneapolis, Minnesota, USA
- 5 Neuroscience Research Institute (NRI) of Gachon University of Medicine and Science, Incheon, South Korea
- 6 Advanced Imaging Research Center (AIRC), Oregon Health & Science University, Portland, Oregon, USA
- 7 Erwin L. Hahn Institute for Magnetic Resonance Imaging (ELH), Essen, Germany
- 8 Center for Imaging in Biomedicine (CIBM), École polytechnique fédérale de Lausanne (EPFL), Lausanne, Switzerland
- 9 Max Planck Institute for Biological Cybernetics (MPI KYB), Tübingen, Germany (9.4T)
- 10 NeuroSpin, French Alternative Energies and Atomic Energy Commission (CEA), Saclay, France
- 11 NeuroSpin, French Alternative Energies and Atomic Energy Commission (CEA), Saclay, France (11.7T)
- 12 Magnetic Resonance Research Center (MRRC), University of Pittsburgh Medical Center (UPMC), Pittsburgh, Pennsylvania, USA
- 13 Max Planck Institute for Human Cognitive and Brain Sciences (MPI), Leipzig, Germany
- 14 Excellence Center for Highfield MR, Medical University of Vienna (MUW), Vienna, Austria
- 15 German Cancer Research Center (DKFZ), Heidelberg, Germany
- 16 Institute of Neuroscience and Medicine (INM), Research Centre Jülich, Jülich, Germany (9.4T)
- 17 Center For Magnetic Resonance And Optical Imaging (MMRCC), University of Pennsylvania Health System (HUP), Philadelphia, Pennsylvania, USA
- 18 Berlin Ultrahigh Field Facility (B.U.F.F.), Experimental and Clinical Research Center (ECRC), Berlin, Germany
- 19 State Key Laboratory of Brain and Cognitive Science, Institute of Biophysics, Chinese Academy of Sciences (CAS), Beijing, China
- 20 Oxford Centre for Functional MRI of the Brain (FMRIB), University of Oxford, Oxford, UK



- 21 Magnetic Resonance Imaging Research Center, Auburn University, Auburn, Alabama, USA
- 22 Center for MR Research (CMRR), University of Minnesota, Minneapolis, Minnesota, USA
- 23 Functional MRI Facility (FMRIF), National Institute of Mental Health and Neurological Disorders and Stroke, National Institutes of Health (NIH-NIMH & NINDS), Bethesda, Maryland, USA
- 24 National Institute of Neurological Disorders and Stroke, National Institutes of Health (NIH-NINDS), Bethesda, Maryland, USA (11.7T)
- 25 National Institute of Information and Communication Technology (NiCT) / Center for Information and Neural Networks (CiNET), Osaka, Japan
- 26 Center for MR Research (CMRR), University of Minnesota, Minneapolis, Minnesota, USA (10.5T)
- 27 Center for Imaging of Neurodegenerative Diseases (CIND), San Francisco VA Medical Center, UCSF, San Francisco, California, USA
- 28 German Center for Neurodegenerative Diseases (DZNE), Bonn, Germany
- 29 Biomedical Research Imaging Center (BRIC), University of North Carolina (UNC), Chapel Hill, North Carolina, USA
- 30 Maastricht Brain Imaging Centre (M-BIC), Maastricht University, Maastricht, The Netherlands
- 31 Maastricht Brain Imaging Centre (M-BIC), Maastricht University, Maastricht, The Netherlands (9.4T)
- 32 Mt Sinai School of Medicine, New York City, New York, USA
- 33 Cleveland Clinic, Cleveland, Ohio, USA
- 34 Centre for Advanced Imaging, University of Queensland, Brisbane, Queensland, Australia
- 35 Royal Melbourne Hospital, University of Melbourne, Victoria, Australia
- 36 University of Sao Paulo (USP), Sao Paulo, Brazil
- 37 Centre d'Exploration Métabolique par Résonance Magnétique (CIMEREM), Marseille, France
- 38 Centre for Functional and Metabolic Mapping, Robarts Research Institute, London, Ontario, Canada
- 39 National Institute for Physiological Sciences (NIPS), Okazaki, Japan
- 40 Kyoto University, Kyoto, Japan
- 41 Zhejiang University, Hangzhou, China
- 42 Brigham and Women's Hospital (BWH), Boston, USA
- 43 University of Southern California (USC), Los Angeles, California, USA
- 44 Cardiff University Brain Research Imaging Centre (CUBRIC) Cardiff, UK
- 45 Wolfson Brain Imaging Centre (WBIC), University of Cambridge, Cambridge, UK
- 46 Imaging Centre of Excellence (ICE), South Glasgow University Hospital, Glasgow, UK
- 47 Magnetic Resonance Research Center (MRRC), Yale University, New Haven, Connecticut, USA
- 48 Comprehensive Heart Failure Center (CHFC), Würzburg University Hospital, Würzburg, Germany
- 49 Weizmann Institute of Science, Rehovot, Israel
- 50 Mayo Clinic, Rochester, USA
- 51 Toronto Western Hospital (TWH), University Health Network (UHN), Toronto, Canada
- 52 National Institute of Health, National Institute on Drug Abuse (NIH-NIDA), Bethesda, Maryland, USA
- 53 Forschungszentrum Jülich, Jülich, Germany
- 54 CRC, University of Liege, Liege, Belgium
- 55 Houston Methodist, Houston, Texas, US
- 56 Athinoula A. Martinos Center for Biomedical Imaging of MGH, Boston, Massachusetts, USA
- 57 University Clinic Erlangen, Erlangen, Germany
- 58 Montreal Neurological Institute and Hospital (MNI), McGill University, Montreal, Canada
- 59 Sungkyunkwan University (SKKU), Seoul, South Korea
- 60 Balgrist Hospital, Zürich, Switzerland
- 61 King's College London (KCL), London, UK
- 62 Barnes-Jewish and Children's Hospital (BJC), St. Louis, Missouri, USA
- 63 Fudan University, Fudan, China
- 64 Swiss Institute for Translational and Entrepreneurial Medicine and Inselspital Bern (item-insel), Bern, Switzerland
- 65 University of California Berkeley, Berkeley, California, USA
- 66 Norwegian University of Science and Technology (NTNU), Trondheim, Norway
- 67 Wellcome Centre for Human Neuroimaging, University College London, UK
- 68 Centre Hospitalier Universitaire de Poitiers (CHU), Poitiers, France
- 69 Zhongnan University Xiangya Hospital, Hunan, China

## Service and exchange – Comprehensive services



Siemens' end-to-end services ensure you stay at the leading edge of MRI technology throughout the entire system lifecycle – from installation, to operation, to upgrades, to ongoing support. Moreover, our diverse communication platforms and communities keep you up to speed on the world of MRI and enable you to share your ideas and experiences with your peers.





#### **Utilization management and reporting**

This powerful solution gives you more from your MRI scanner. It allows you to monitor KPIs and benchmark your system against other Siemens MRI machines at any facility or organization. So you can keep track of your MRI performance, and reap the maximum reward from your scanner.

#### **Predictive maintenance**

When systems go down, it impacts both your ability to care for your patients and your bottom line. Siemens provides a predictive maintenance service to help you minimize lost time. It informs you when a part of your MRI system is likely to fail, enabling you to plan repairs and prevent downtime before it happens.

#### **EVOLVE**

Keep your hardware and software up to date at all times – a key factor in enhancing performance and diagnostic quality. You receive all applicable upgrades for software and the *syngo* OS, plus at least one workstation hardware upgrade within the first six years.


#### **Siemens Guardian program**

This program provides the latest service technology so you can better manage your MRI system. It combines many features in a single package – offering real-time system monitoring, expert advice to improve workflow efficiency, proactive maintenance, and support. Moreover, it guarantees defined repair times, giving you complete peace of mind.

#### **Proven upgrade paths**

With MAGNETOM scanners, taking your MRI system to the next level is simplicity itself, thanks to clearly defined upgrade paths. In fact, Siemens has built an entire organization (CDV) to help customers truly maximize their system life – and increase their return on investment.

## Service and exchange – Peer-to-peer information



*On MAGNETOM Flash:*

*“An excellent and useful combination of technological and clinical articles that both keep one up to date with advances in MRI and provide practical assistance for day-to-day practice – good and interesting learning material.”<sup>13</sup>*

Mark Lourensz  
St Vincent's Hospital, Fitzroy, Victoria, Australia



#### **MAGNETOM World**

Siemens Healthineers' global MRI community offers peer-to-peer support and information. Radiologists, cardiologists, technologists, and physicists have all contributed with publications, presentations, training documents, case studies, and more – all freely available to you via this unique network. Plus, the bi-annual MAGNETOM World Summit is the ideal opportunity to share and exchange ideas.

#### **MAGNETOM Flash**

Published quarterly, the MR customer magazine features up-to-date clinical case studies, application tips and technical and product information relevant to you. All content is carefully compiled by experts to meet the needs of today's MRI users in both clinical and research scenarios. In fact, 98.5% of readers report that MAGNETOM Flash is clinically relevant.

#### **IDEA / UHF Online Discussion Board**

IDEA<sup>8</sup> is an open development platform for the largest and most active 3T and UHF research communities in the world. It unites users from across the globe and fosters innovation in the field of MRI. Members collaborate online at [www.mr-idea.com](http://www.mr-idea.com) and at an annual meeting. IDEA includes an exclusive area, the UHF Online Discussion Board, to help users focus on topics of interest, as well as find and communicate with the right peers.

**Visit MAGNETOM World**

[siemens.com/magnetom-world](http://siemens.com/magnetom-world)

# Technical specifications

## MAGNETOM Terra Technical specifications

Field strength	7 Tesla
Bore size	60 cm
Magnet length	270 cm
System length	297 cm
System weight (in operation)	< 25 tons
Minimum room size <sup>17</sup>	65 m <sup>2</sup> / (w/o pTX <sup>2</sup> and 3rd order shim option)
Dual Mode functionality <sup>2</sup>	Clinical Mode Research Mode <sup>2</sup>
RF transmit	TimTx-1, TimTx-8, TimTx-16 (in research mode) <sup>2</sup>
Maximum number of channels <sup>3</sup>	32, 64
Number of independent receiver channels that can be used simultaneously in one single scan and in one single FOV, each generating an independent partial image	32, 64
Gradient strength	XR gradients (80 mT/m @ 200 T/m/s)
Helium consumption	Zero Helium boil-off technology
Quiet Suite/DotGO	Available in research mode <sup>2</sup>
Local coils in clinical mode	1TX/32RX head coil 1TX/28RX knee coil 23Na 1Tx/32Rx head coil 31P/1H TxRx flex loop coil





## Why Siemens Healthineers?

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey towards expanding precision medicine, transforming care delivery, and improving patient experience, all enabled by digitalizing healthcare.

An estimated 5 million patients globally everyday benefit from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics and molecular medicine, as well as digital health and enterprise services<sup>18</sup>.

We are a leading medical technology company with over 170 years of experience and 18,000 patents globally. With more than 48,000 dedicated colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.



Not for distribution in the USA. On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States. Some products are still under development and not commercially available yet. Their future availability cannot be ensured. The information in this document contains general technical descriptions of specifications and optional features which do not always have to be present in individual cases. Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

For accessories, please visit:  
[www.siemens.com/medical-accessories](http://www.siemens.com/medical-accessories)

- 1 Compared to 3T systems.
- 2 Research mode as part of dual mode is available as an option and not intended for clinical use. Research operation may require observation of national regulations.
- 3 Channels (coil elements) that can be connected simultaneously.
- 4 Keil et al., Magn Reson Med 70:248–258 (2013); Wiesinger et al., Magn Res Med 52:953–964 (2004); Pruessmann et al., Magn Reson Med 42:952–962 (1999); Griswold et al., Magn Reson Med 47:1202–1210 (2002)
- 5 Compared to previous 7T generation.
- 6 Under normal operating conditions with standard Siemens sequences/protocols.
- 7 Example images available in this Brochure: Page 12; SWI minIP/phase, Page 13; T2 TSE, SWI, Page 15, Page 16; PD FSE TSE, Page 17; SWI, T2 TSE, Page 22 ; SWI, Page 23 ; PD TSE FS, Page 24; T1 SE, Page 25 ; PD qTSE FS.
- 8 Scheenen et al., Magn Reson Mater Phy 21:95–101 (2008)
- 9 Heidemann et al., Magn Reson Med 68:1506–1516 (2012); Yacoub et al., PNAS 105:10607–10612 (2008)
- 10 Madelin et al., J Magn Reson Imaging 38:511–529 (2013); Valkovic et al., Analytical Biochemistry 529:193–215 (2017)
- 11 <https://health.usnews.com/best-hospitals>
- 12 <https://www.ismrm.org/18m>
- 13 The statements by Siemens' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.
- 14 Trattnig, et al., NMR Biomed. 9:1316–1334 (2015)
- 15 Cinematic VRT is recommended for communication, education, and publication purposes and not intended for diagnostic reading.
- 16 Rendered with a Siemens internal cinematic rendering prototype.
- 17 Minimum total space requirement for magnet, electronics, and console room.
- 18 Siemens AG, "Sustainable healthcare strategy – Indicators in fiscal 2014", pages 3–4.

.....

## Siemens Healthineers Headquarters

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
[siemens-healthineers.com](http://siemens-healthineers.com)